



# Experience the Trace Quality



*Happier Living Everyday*

## SINGLE CHANNEL ECG RECORDER **CARDIART 6108T**



### Single Channel ECG Recording

Enables the recording of all leads sequentially



### Ergonomic Design

The compact, light weight and portable design ensures ease of use



### Automatic and Manual Recording Modes

Manual mode enables the acquisition and printing of ECG recording in real time. Simultaneous acquisition is possible with automatic mode



### Last Memorised ECG Recording

Ensures the storage of last recording in the memory



### Thermal Printing

Allows the printing of ECG on 50 mm paper using a high-resolution thermal print



### Power Supply

Comes with dual power supply – mains and rechargeable internal battery



### Built-in Rechargeable Battery

Enables mains independent operations and provides battery use after charging

#### CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited  
11th KM, Bannerghatta Road,  
Arakere, Bangalore - 560076  
CIN: U33110KA2012PTC067282  
Ph: +91 80 26484388/ 2648 4350  
Toll Free: 1800-4252355  
Website: [www.bplmedicaltechnologies.com](http://www.bplmedicaltechnologies.com)  
For Enquiries: [sales.medical@bpl.in](mailto:sales.medical@bpl.in)



A BPL Medical Technologies  
Magazine



Follow Us on

