

With
APAP



CPAP & APAP

BPL *Harmony 1*



Happier Living Everyday



Multiple therapy modes

Comes with CPAP and APAP therapy modes with auto ramp function



Humidifier

Equipped with an in-built humidifier along with a pre-heating function



Intelligent Pressure release (iPR)

Makes the therapy comfortable, by reducing the treatment pressure upon exhalation



Power efficient

Equipped with energy saving mode and an auto power on/off thereby improving the machine durability



Smart Solutions

Has high-resolution SD card data storage
Wireless transmission enabled through Wi-Fi, Bluetooth and GPRS*



Intelligent adjustments

BPL Harmony comes with auto altitude adjustment and auto screen luminance adjustment*



Follow Us on



*optional

CPAP MACHINE

BPL HARMONY

Product Specifications

Display	3.5" Colour LCD DISPLAY
Flow	>120 L/min
Safety specifications	Class II Type BF
Ingress Protection	IP22
Air filter Efficiency	>75% (7-micron dust)
Pressure Setting Range Max Single Fault Steady Pressure Pressure Control Accuracy Dynamic Pressure Variation Pressure monitoring Accuracy	4-20 cmH2O 40 cmH2O ± 0.5 cmH2O < 1 cmH2O $\pm (2\% \text{ full scale reading} + 4\% \text{ actual reading})$
Humidifier Water capacity Humidity	290 ml (MAX Water Level) >10 mg/L BTPS (Within the set pressure range)
Noise Value A-weighted sound pressure level A-weighted sound power level	≤ 28 dBA ≤ 36 dBA
Dimension Dimension (L*W*H) Weight	235*148*125 mm Approximately 1.35 kg
Electrical Specifications AC Power Input DC Power Input	100~240 V~ 50/60 Hz 2.0 A Max 24 V, ≤ 2.5 A
Operating environment conditions Temperature Relative Humidity Atmospheric Pressure	5°C~35 °C 10%~95%(non-condensing) 70 kPa~106 kPa
Storage Environment conditions Temperature Relative Humidity Atmospheric Pressure	-25°C~70°C 5%~95% (non-condensing) 70 kPa~106 kPa

CERTIFIED ISO 13485 : 2016 COMPANY

BPL Medical Technologies Private Limited
 Regd. Office: 11th KM, Bannerghatta Road,
 Arakere, Bangalore - 560076, India.
 Toll Free: 1800-4252355
 Website: www.bplmedicaltechnologies.com
 For Enquiries: sales.medical@bpl.in
 CIN: U33110KA2012PTC067282

