



Experience the Finest Display

MULTI PARAMETER MONITOR
ULTIMA PRIME



Happier Living Everyday



**Provides connectivity to
BPL Medius+ Central Nursing
Station Software and displays a
maximum of 10 waveforms**



High Resolution Display

High resolution 12.1" TFT colour display,
optional Touch screen



Standard Configuration

Standard configuration of ECG, RESP, Nellcor SpO₂ /
Masimo SpO₂ / Digital SpO₂, NIBP, TEMP



Calculation

The added functions of the device are Drug dose
calculation, HRV analysis & Hemodynamic calculation.



Connectivity to CNS

Connectivity to optional BPL Medius+ Central
Nursing Station Software (Wired or optional wireless)



168 Hour Trend

Provides for maximum of 168 hours of trend review



ECG Monitoring

5 Lead ECG / 3 Lead ECG (optional)



Waveform Display

Maximum 10 waveforms display



Arrhythmia and ST Analysis

The device provides for arrhythmia detection and
ST analysis



Optional Parameters & Configuration

2-IBP, ETCO₂ (Main Stream / Side Stream),
Multi-Gas / AGM, Noninvasive Cardiac Output, HL7
Output and Thermal Recorder