





RECORDERS & MEDICARE SYSTEMS

Plot # 196, Industrial Area, Phase-1, Panchkula, Haryana INDIA - 134113

Patient: UNM
Refd. By: Dr. R.L.
Pred. Eqns: RECORDERS
Date : 25-Dec-1999 11:25 AM

Age : 56 Yrs 04 Mths Gender : Male
Height : 167 Cms Smoker : No
Weight : 62 Kgs Eth. Corr: 100
ID : 1 Temp : 15°C



Spirometry Results

Parameter	Pred	M.Pre	%Pred	M.Post	%Pred	%Imp
FVC (L)	03.19	03.24	102	03.44	108	+06
FEV1 (L)	02.54	02.84	112	02.85	112	---
FEV1/FVC (%)	79.62	87.65	110	82.85	104	-05
FEF25-75 (L/s)	03.52	02.97	084	02.78	079	-06
PEFR (L/s)	08.36	07.43	089	07.38	088	-01
FIVC (L)	---	03.24	---	03.39	---	+05
FEV.5 (L)	---	02.15	---	02.21	---	+03
FEV3 (L)	03.09	03.24	105	03.44	111	+06
PIFR (L/s)	---	06.64	---	06.94	---	+05
FEF75-85 (L/s)	---	01.22	---	00.88	---	-28
FEF.2-1.2 (L/s)	06.19	06.16	100	06.52	105	+06
FEF 25% (L/s)	07.62	06.16	081	06.69	088	+09
FEF 50% (L/s)	05.30	03.31	062	03.08	058	-07
FEF 75% (L/s)	02.35	01.60	068	01.23	052	-23
FEV.5/FVC (%)	---	66.36	---	64.24	---	-03
FEV3/FVC (%)	96.87	100.00	103	100.00	103	---
FET (Sec)	---	01.79	---	02.97	---	---
ExptTime (Sec)	---	00.05	---	00.05	---	---
Lung Age (Yrs)	048	042	088	042	088	---
FEV6 (L)	03.19	---	---	---	---	---
FIF 25% (L/s)	---	05.92	---	06.75	---	+14
FIF 50% (L/s)	---	06.54	---	06.49	---	-01
FIF 75% (L/s)	---	05.48	---	04.59	---	-16
SVC (L)	---	03.49	---	03.13	---	-10
ERV (L)	01.24	00.85	069	00.78	063	-08
IRV (L)	---	01.76	---	01.42	---	-19
VE (L/min)	---	22.96	---	26.57	---	+16
Rf (l/min)	---	26.09	---	28.57	---	+10
Ti (sec)	---	01.00	---	00.90	---	-10
Te (sec)	---	01.30	---	01.20	---	-08
VT (L)	---	00.88	---	00.93	---	+06
VT/Ti	---	00.88	---	01.03	---	+17
Ti/Ttot	---	00.43	---	00.43	---	---
IC (L)	---	02.64	---	02.35	---	-11
MVV (L/min)	123	102	083	121	098	+19
MRf (l/min)	---	54.55	---	68.57	---	+26
MVT (L)	---	01.86	---	01.76	---	-05

Pre Medication Report Indicates

Spirometry within normal limits as (FEV1/FVC)%Pred >95 and FVC%Pred >80

Post Medication Report Indicates

Spirometry within normal limits as (FEV1/FVC)%Pred >95 and FVC%Pred >80

Doctor's Notes

Patient has responded well

Dr. M.K