



Experience the Trace Quality



SIX CHANNEL ECG RECORDER CARDIART 8108 VIEW

Happier Living Everyday



Six Channel ECG Recording

Simultaneous 6 channel ECG recording with 12 lead simultaneous acquisition



LCD Display

3.5 inch LCD display with One-Touch operation enables the easy viewing of patient demographics and real time ECG data option



Alphanumeric Keyboard

Ensures the freedom of ease of use with one hand



Easy Measurement and Interpretation

Built-in ECG parameters measurement and interpretation



Built-in Rechargeable Battery

High performance Li-ion rechargeable battery enables mains independent operations and provides convenience of easy use



Print Mode

Different print modes include pre-sample, real-time sample, period sample and arrhythmia triggered sample



Internal Memory

Has an Internal Memory, which can store up to 800 ECG recordings



PC Interface

Interfaces with PC using optional Smart-View software