



Experience the Quality Care

FOETAL MONITOR

DIVA - FM 9852/ 9853/ 9854



Happier Living Everyday



Offers clear and high resolution display that helps to make more accurate decisions

FM 9852



Monitor Screen

7" high resolution colour TFT display with tiltable screen



Standard Configuration

FHR, TOCO, and foetal movement



Portability

The device is compact, portable and lightweight



Display Interface

Large letters display interface helps in making accurate decisions even from a distance



Ultrasound Probe

The probe with high sensitivity enables user to monitor foetus > 12 weeks and provide versatility in clinical settings



Data Storage

The device can record data for 12 hours and comes with playback and print facilities



Thermal Printer

Thermal printer facilitates easy availability of printed results



Safety

The low frequency ultrasound waves are safe on foetus



AFM

Automatic Foetal Movement detection with event marker



Li-ion Battery

The built-in rechargeable Li-ion battery ensures uninterrupted use of the device



Central Monitoring System Software

This is an optional feature.

*Variant

Twin FHR monitoring facility is available with FM 9852T



Provides an audible simulation of the foetal bloodflow in more accurate and sensitive manner

FM 9853/ 9854



Monitor Screen

12.1" high resolution colour TFT display with tiltable screen upto 90°



Standard Configuration

FHR, TOCO, and foetal movement



Portability

The device is compact, portable and lightweight



Ultrasound Probe

The probe with high sensitivity enables user to monitors foetus >12 weeks and provide versatility in clinical settings



Data Display & Storage

The device can record and display data for 25 hours with dynamic data save capability



Thermal Printer

Thermal printer facilitates easy availability of printed results



Safety

The low frequency ultrasound waves are safe on foetus



AFM

Automatic Foetal Movement detection with event marker



Li-ion Battery

The built-in rechargeable Li-ion battery ensures uninterrupted use of the device



Central Monitoring System Software

This is an optional feature.

*Variant

Twin FHR Monitoring is available in FM 9854

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies Magazine



Follow Us on

