



## D100 Nasal CPAP Mask

### Simple, Adjustable Comfort

The D100 Nasal CPAP Mask offers comfort, durability and simplicity for effective CPAP therapy while remaining cost-effective for both patients and providers alike. Available in three sizes, the D100 Nasal mask is extremely easy to fit, clean and maintain. It can be easily adjusted using the headgear straps to achieve an effective seal, leading to improved patient compliance. The D100 Mask is also available as a full face mask.



D100 Nasal CPAP Mask

### Features

- Four-point adjustable headgear for a secure, customized fit
- Lightweight, molded silicone cushion allows for comfortable extended use
- Three sizes of forehead pad included with every mask for maximum comfort
- 360° swivel for optimal tube positioning and easier mobility

## Ordering Information

### Product Description

#### Complete Mask

Small  
Medium  
Large

### Item Number

D100N-S  
D100N-M  
D100N-L

#### Replacement Mask Cushions

Small  
Medium  
Large

D100N-CS  
D100N-CM  
D100N-CL

### Product Description

#### Replacement Silicone Forehead Pads

Small  
Medium  
Large

### Item Number

D100-PS  
D100-PM  
D100-PL

#### Parts & Accessories

Headgear

D100-H



The IntelliPAP® System is the perfect fit with the D100 Nasal Mask.

## Reimbursement Information

Product Description	Code Description	HCPSC Code	Frequency	Allowable Range
D100 Nasal Mask	Nasal Application Device (mask or cannula type), used with positive airway pressure device, with or without head strap	A7034	1 per 3 months	\$109.77-\$131.75
Replacement Headgear	Headgear used with positive airway pressure device	A7035	1 per 6 months	\$32.37-\$38.85

The noted allowable rates are based off of the January 2013 HCPSC Fee Schedule. All information presented is for reference purposes only.



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 | USA  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)