

565 Nucla Way, Unit B Aurora, Colorado 80011 USA

## Nuprep Gel Skin Prep Gel





## **About Nuprep Skin Prep Gel**

For quality tracings, proper skin preparation at every electrode site is crucial.

Nuprep Skin Prep Gel gently abrades and moistens the top layer of skin. This lowers the skin impedance with minimal skin irritation and patient discomfort.

Nuprep virtually eliminates the undesirable effects of diaphoresis and unwanted artifact, thus providing the best opportunity to achieve accurate neurodiagnostic and cardiodiagnostic tracings.

## **Target Markets**

- **EEG Examinations**
- Sleep Diagnostics
- Evoked Potential Procedures
  - Intraoperative Monitoring
    - Audiology
    - Biofeedback
    - Cardiac Stress Tests
      - Holter Monitoring
- Cardiac Rehabilitation Monitoring

## **Contact Us**

(800)525-2130

www.weaverandcompany.com info@weaverandcompany.com