



Positive Airway Pressure Device

Easy & Effective

Comfortable | Sleep

Respirare Technologies



Integrated Humidification

Heated Humidification makes therapy more comfortable and can reduce side effects like a dry nose and throat.



Ex.t

The Ex.t is a comfort feature that will lower the pressure when you exhale and slowly increase the pressure when you inhale and will be available in APAP mode.



Key Features

Supported by RespiView

RespiView is Deck Mount's PC-based clinical analysis and patient data management software. It enables you to review device therapy settings and download, analyze and store therapy data.

Ramp for extra comfort

Respirare CPAP offers Ramp feature to help patients ease into each therapy session. Ramp reaches the prescribed therapy level gradually and comfortably.

Product Codes

Respirare APAP with Humidifier

VT 40

Note: Technical specification may change without notice. Always refer to your device user manual.

Technical Specification

Modes	CPAP , APAP	
Operating Range	4-20 CmH ² O	
Sound Pressure Level	<28 dBA with uncertainty of 2 dBA	
Leak Compensation	10-40 LPM	
Peak Flow	200 LPM	
Ramp Time	Auto, 0-60 mins	
Ex.t (Expiratory Pressure relief)	1, 2, 3	
Humidification Level	Off, 1-5 level	
Alarms	Mask Alert	
Data Management	Micro SD card	
Dimensions (LxWxH)	290x160x110 mm (Device with Humidifier)	
Weight	1800 g (Device with Humidifier)	
Power Adapter	Input Range : 150-240V, 50/60 H ² , 5A 60 Watt (typical Power Consumption)	
Operating Temperature/ Humidity/ Altitude	+5 c to +35 c/ 10-95% relative humidity, non-condensing/sea level to 2,591 m, air pressure range 1013 hpa to 735 hpa	
Storage and Transport Temperature/Humidity	-20 c to +60 c / 5-95% relative humidity, non-condensing	
Supplement Oxygen	Maximum Flow: 4L/min (CPAP)	
Standard Air Filter	Pollen Filter (Washable)	
Tubing	Standard (19 mm)	
Air Outlet	22 mm Air Outlet (STD)	
Warranty	2 Years (Replacement)*	
IEC	IEC 60601-1-1 IEC 60601-1-2 IEC 80601-2-12	Class 2 (double insulated, Type BF, Ingress Protection IP22)

*Term & Conditions Apply

Respirare Technologies



Scan Report

Simplest way to get one day compliance report by scanning the QR code.



Gender Selection

Gender Selection helps to give a Soft/Standard response to the selected patient accordingly.

Technical Specification

Modes	CPAP , Auto CPAP	
Operating Range	4-20 CmH ² O	
Sound Pressure Level	<28 dBA with uncertainty of 2 dBA	
Leak Compensation	10-40 LPM	
Peak Flow	200 LPM	
Ramp Time	Auto, 0-60 mins	
A-Flex (Expiratory Pressure relief)	1, 2, 3	
Gender Selection	Male/ Female	
Humidification Level	Off, 1-5 level	
Alarms	Mask & Leak Alert	
Data Management	Yes (1 day, 7 day, 31 days) of detailed data stored in the device Scan to see one day compliance report	
Dimensions (LxWxH)	290x160x110 mm (Device with Humidifier)	
Weight	1800 g (Device with Humidifier)	
Power Adapter	Input Range : 150-240V, 50/60 H ² , 5A 60 Watt (typical Power Consumption)	
Operating Temperature/ Humidity/ Altitude	+5 c to +35 c / 10-95% relative humidity, non-condensing/sea level to 2,591 m, air pressure range 1013 hpa to 735 hpa	
Storage and Transport Temperature/Humidity	-20 c to +60 c / 5-95% relative humidity, non-condensing	
Supplement Oxygen	Maximum Flow: 4L/min (CPAP)	
Standard Air Filter	Pollen Filter (Washable)	
Tubing	Standard (19 mm)	
Air Outlet	22 mm Air Outlet (STD)	
Warranty	2 Years (Replacement)*	
IEC	IEC 60601-1-1 IEC 60601-1-2 IEC 80601-2-12	Class 2 (double insulated, Type BF, Ingress Protection IP22)

*Term & Conditions Apply

Key Features

A-Flex

The A-Flex is a comfort feature that will lower the pressure when you exhale and slowly increase the pressure when you inhale and will be available in Auto CPAP mode.

i-Mode

This feature will turn the device **On** for you have the mask on without the need to press any buttons and will also turn the device **Off** once you take your mask Off.

Integrated Humidification

Heated Humidification makes therapy more comfortable and can reduce side effects like a dry nose and throat.

Product Codes

Respirare Auto CPAP with Humidifier

VT 50

Note: Technical specification may change without notice. Always refer to your device user manual.

Manufactured By:
Deck Mount Electronics Pvt. Ltd 260, Phase IV,
Udyog Vihar, Sector 18, Gurugram, Haryana 122015

Saving millions of lives

DeckMount Respirare CPAP with A-Flex offers the ultimate sleep therapy comfort. It combines auto-adjusting pressure delivery with breath-by-breath pressure relief for most natural sleep.



wecare@deckmount.in



www.deckmount.in



Toll Free No.
1800 30 92 499



Deck Mount Electronics
Pvt. Ltd, 260, Udyog Vihar,
Phase 4, Gurugram 122015,
Haryana, India



Make In India

