



BPL *LifePAP* series
BPL *Harmony* series



Happier Living Everyday



Follow us on
    

www.bplmedicaltechnologies.com

BPL *LifePAP* series

BPL *Harmony* series



Stable quality equipped with highly advanced algorithm



High Performance Motor

Silent DC brushless motor, 0.2 seconds quick response, 20000 hours of service life.



No spill Air Duct design

Accurately monitor the air flow signal to give seamless experience to the user.



Highly Accurate sensors


Equipped with advanced pressure and flow sensors, Accurately collect pressure & flow signal changes and provide real-time control to users.




Multiple Smart Alerts

Equipped with smart alerts such as High Pressure alert, Leak Reminders etc to provide a safe and convenient use.

Excellent User & Device Synchrony



Fast & Accurate response | Optimal treatment effects | High sensitivity



Smart Detections
If the air leakage reaches 60lpm, it can still be triggered synchronously.

User Friendly Design


01 Easy Separated Structure

Humidifier can be separated from the device, the host can be used individually.



03 Data management & Smart Analysis

Equipped with PC APP for patient data management and smart analysis



02 Convenient Operated Button

The simple combination of two buttons finish all operations.



04 Lock Type

Ergonomic design constructed with Medical grade PC material. Safe and secure to use and operate



Specification	BreathCare PAP	
Picture	 <p>White&Gray</p>	
Model	Harmony Prime	LifePAP PRIME
Mode	CPAP, Auto CPAP	CPAP, S, T, ST, VGPS
Pressure	4-20cmH ₂ O	4-25cmH ₂ O
Efficient humidifier	●	●
Smart Start/Stop	●	×
FPS-Tech	●	×
High Pressure Reminder	×	●
Leak Reminder	●	●
Asphyxiation Reminder	×	●
Trigger Adjustment	×	●
Cycle Adjustment	×	●
Slope Adjustment	×	●
Tidal Volume Adjustment	×	●
Altitude Adjustment	●	●
Leak Compensation	●	●
Data Management	●	●

ISO 13485:2016 CERTIFIED COMPANY
ISO 9001:2015 CERTIFIED COMPANY

BPL MEDICAL TECHNOLOGIES PRIVATE LIMITED
Regd. Office: 11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076, India.
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in
CIN: U33110KA2012PTC067282



Happier Living Everyday

