

POWER BACKUP FOR CPAP & BIPAP





POWER BACKUP FOR CPAP & BIPAP

Power backup for CPAPA/BIPAP machines, compatible with multiple brands like resmed, philips, respicare, hcl etc.

USB-C 60W | DC 12-24V output | 10W Wireless Fast Charging



COMES IN 2 DIFFERENT OPTIONS



(With 4 Hours Of Back Up
For CPAP&BIPAP)



(With 8 Hours Of Back Up
For CPAP&BIPAP)

COMPATIBLE CPAP MODELS

(Not Limited to the models below)



Phillip Dreamstation



Phillip Respiroics REMstar



Resmed S9



Resmed Airsense 10



COMES WITH 3 DIFFERENT TYPE OF CONNECTORS



- There are three different types of DC connectors, which you can select as per your device input port
- Direct DC cable increases the efficiency of the product by reducing the hassle of connecting adapters

MULTIPLE CONNECTORS



CPAP MACHINE



BIPAP MACHINE



BIPAP MACHINE



SPECIFICATIONS

Backup	4 Hrs & 8 Hrs
Model No.	CMPBXIP-56
Battery Type	Lithium-ion
Input	DC : 18-24V / 3A USB-C : 5V/3A, 9V/3A, 12V/3A, 15V/3A, 20V/3A
Output	USB-C : 5V/3A, 9V/3A, 12V/3A, 15V/3A, 20V/3A USB-C2 : 5V/3A, 9V/3A, 12V/2.25A USB-A : 4.5V/5A, 5V/4.5A, 9V/2A, 12V/1.5A 12V/4A, 16.5V/4A, 20V/4A, 24V/3.5A
DC Output	12V/4A, 16.5V/4A, 20V/4A, 24V/3.5A
Wireless Output	10W
Working Temperature	0°C ~ 60°C
Display	Power Indicated by LCD
Charge Retention	>5000h
Material	ABS+PC
Charging Hours	1-3 Hours
Warranty	5 Years
Size	122*122*26.5 mm
Weight	500 gms.