

**LUCKCOME**

Care Life From The Beginning

## L8 Fetal/Maternal Monitor

7" TFT display



L8 TFT Display



L8 LED Display



L8 LCD+LED Display

## L8 FETAL /MATERNAL MONITOR

### 7" TFT display

#### Technical Specifications

#### FEATURES:

- 7" TFT color screen ,high resolution. Folding up 90 degree
- Build in long life 110mm thermal printer
- Special high sensitive probe .
- Long time data storage which can be played back and printed
- Rechargeable battery many hours work
- DSP technology. Real time identifying and measurement for FHR ,accurate and reliable results
- Automatic Fetal Movement Detection and Accounting
- Central Monitoring System
- CTG Scoring sytern:Krebs score/Fischer Score

#### SPECIFICATION:

##### FHR

Transducer: Multi-crystals, wide wave beam,pulsed-wave working mode

Strength: < 5mW/cm<sup>2</sup>

Working Frequency: 1.0MHz

##### TOCO

Measurement Range: 0-100 units

Display: 7" TFT LCD

Dimension: 295D×240W×73H(mm)

Net Weight: 1.75kg

#### OPTIONS:

- Touchscreen function
- Connect with L8T Wireless Station
- L8A (Maternal SpO<sub>2</sub>, PR, NIBP)
- FAS (Fetal Acoustic Stimulator)
- Twins FHR Measurement