

OXYMED[®]

New RE Series

Experience market leading PAP devices-
Empowering precise clinical decisions with
real time SpO₂ monitoring



Safety



Comfort



Ultra-quiet
Operation



Large Vivid
LCD Display



PRODUCT INTRODUCTION



Introducing New generation PAP devices, with intelligent data management and enhanced comfort redefined.

EFFICIENCY



Wide pressure adjustment range ensures optimal therapy for different users

Integrated data management system supports clear monitoring of therapy progress

COMFORT



Climate control- Intelligent temperature and humidity control system

AST & COMF technology to ensure optimal comfort and therapy efficiency

SAFETY



Bluetooth SpO₂ monitoring to ensure users better clinical outcomes.

Automatic detecting 6 events and responding to 4 events to ensure optimal treatment pressure

DESIGN




Ergonomically designed 3.5 inches tilted display for enhanced visibility and an improved user experience

User data information showed on the screen

TECHNICAL SPECIFICATION

Model	SleepEasy RE	Airsart Auto RE	Airsart ST25 RE	Airsart ST30 RE
Mode	CPAP AUTO CPAP - -	CPAP AUTO CPAP BPAP AUTO BPAP	CPAP	CPAP, AUTO CPAP, BPAP, AUTO BPAP, S, T, S/T, VAT, PC
Pressure Range	4-20cm H ₂ O	4-25cm H ₂ O	4-25cm H ₂ O	4-30cm H ₂ O
Ramp Time	0 min-45min (5-minute incrementals), Ramp Pressure : 4cmH ₂ O~ (EPAP-0.5) cmH ₂ O			
COMF Pressure Relief	1-3 levels			
Humidification Level	1–5 levels (95 to 167°F / 35 to 75°C)			
Rise Time	-	1-4 levels	1-3 levels	
AHI Display	YES			
Night Staging	YES			
Preheat	5 minutes from 23°C to 45°C			
Auto Drying	YES			
Trigger Sensitivity (S, S/T,PC)	-			Auto, 1-3 levels
Shifting Sensitivity (S,S/T)	-			Auto, 1-3 levels
Trigger Locking (S,S/T)	-			off, 0.3S-1.5S
Blue-tooth	Optional	Standard	Standard	Standard
SpO ₂	Optional	Standard	Standard	Standard
Wi-Fi	Optional			
4G	Optional			

Multiple Data Management Platform



INTELLIGENT DATA MANAGEMENT



REMaster App



Wi-Fi Connection



Routinely check therapy report

REMaster Cloud Platform



User's therapy data is accessible on ResMaster cloud platform



Device and consumables status information are available on ResMaster cloud platform

ACCESSORIES



SD Card



Air filter x3



Full Face Mask



Air Tubing



T-Junction



Power Supply



User Manual



Carrying Bag

OXYMED RE SERIES



MEDEQUIP®
Healthcare Solutions Pvt. Ltd.
Making you Independent

Head Office:

Manufacturing Unit 1 :

Plot No.A-12,A-13, Apparel Park Phase 1, Doddaballapur Industrial Area, Bengaluru-561203, Karnataka, India

Manufacturing Unit 2 :

Plot No A-49 to A-69, Radha Govind Vitaika 3rd,Gram - Kherwadi, Tehsil - Amer, Jaipur - 303805



+91 80 2836 2444 / 2836 2555



customersupport@medequip.co.in



www.medequip.co.in